

**DiaDent**<sup>®</sup>

**DiaDent Group International Inc.**  
[www.diadent.com](http://www.diadent.com)

**The name you trust**

**The product you know**



**Diapex**<sup>®</sup>

Silicone Oil Based Calcium Hydroxide + Iodoform Paste  
for Temporary or Permanent Root Canal Treatment

SAVES TIME

MEDICATED

RADIOPAQUE

MADE IN USA

**DiaDent**<sup>®</sup>

**Corporate Office**

#11 - 3871 North Fraser Way, Burnaby, BC CANADA V5J 5G6

t: 604.451.8851 | f: 604.451.8865 | tf: 1.877.342.3368 | e: [diadent@diadent.com](mailto:diadent@diadent.com) | [www.diadent.com](http://www.diadent.com)

**Easy techniques - Saves time and money**

Eliminate messy time-consuming mixing! Silicone oil based calcium hydroxide paste is pre-mixed with iodoform in a convenient pre-filled syringe. Tips are pre-bent to ensure excellent canal accessibility. This ensures shorter chair time and increases patient comfort.

**Quality of ingredients**

Made in USA, Diapex is made from the purest grade iodoform, increasing Diapex's bacteriostatic and radiopacity properties. Silicone oil ensures complete coating of canal walls and it never hardens so that calcium hydroxide remains active in the root canal, stimulating apexification and apexogenesis. It also allows for ease of manipulation and has superior flow characteristics.



**Diapex**<sup>®</sup> is recommended for the following applications:

- Temporary or permanent root canal treatment
- Vital pulpotomies in deciduous teeth
- Apexification and Apexogenesis
- Stimulation of apexification and hard tissue formation
- Traumatic injuries
- Root resorption treatment including primary teeth

**Techniques for using Diapex**

Easy as 1.2.3...

**Temporary Method:**



1. Clean and shape the canal.
2. Insert the pre-bent tip into the apical fifth of root canal length.
3. Fill the canal.

**Diapex** does not set and can be removed during a subsequent visit using sodium hypochlorite.

**Permanent Method:**



1. Place a small amount of **Diapex** at the periapical portion of the canal to ensure apical closure.
2. Complete obturation using normal Gutta Percha and Sealer techniques.
3. Place a permanent restoration to seal the access opening.

Ask for **Diapex** by name from your local dental supplier.

**Ingredients:**

Calcium Hydroxide	30.3%
Iodoform (purest grade)	40.4%
Silicone Oil (Medical Grade)	22.4%
Inert	6.9%

<b>Regular Kit</b>	<b>(#1001-301):</b>	2.0g of paste + 20 disposable tips
<b>Intro Kit</b>	<b>(#1001-304):</b>	0.5g of paste + 4 disposable tips
<b>Syringe</b>	<b>(#1001-302):</b>	2.0g of paste
<b>Tips</b>	<b>(#1001-303):</b>	20 disposable pre-bent tips
<b>Tips</b>	<b>(#1001-316):</b>	20 disposable straight tips (blister pack) (Thinner tips and can be bent)

**Shelf Life:** 4 years

Diapex was rated as



**Corporate Office**