



The Prelude to your Symphony



Just like a *prelude* to a symphony, Dia-Prep Plus prepares you for the success of your root canal treatment. Dia-Prep Plus is an EDTA cream with Urea Peroxide, used as an irrigation agent for chemo-mechanical preparation of the root canal.

- ♪ Removes smear layer, softens dentin, and facilitates the removal of calcific obstruction
- ♪ Produces oxygen bubbles that flush out debris within the canal
- ♪ Eliminates instrument binding or breakage
- ♪ Permits better penetration of root canal medicaments and fillings



Dia-PrepTM *Plus*

EDTA Root Canal Preparation Cream

DiaDent[®]

DiaDent Group International
www.diadent.com | 1.877.DIA.DENT

Contact your local dealer for product
and pricing information