

DiaDent[®]

1.2 Mega Pixel USB 2.0 Intraoral Camera

DiaCam



Excellent resolution

DiaCam displays 1280x960 (1.2 megapixels) high-resolution digital images.

Easy installation

The USB cable goes into the camera hand-piece and into any USB 2.0 port of a computer after setting up the USB driver and application software (option).

Effective power control

Pushing the power on/off button on the hand-piece enables the camera to turn on/off power or LED. No separate power supply is required. Power will automatically shut down after 5 minutes when the camera is not used.

Complete focus control

DiaCam has the focusing button on the hand-piece to adjust the focus of the lens.

Convenient capture function

DiaCam also has the capture button on the hand-piece to capture images. The captured images can be saved onto computer through the application software.

Visit us at the
ADA Show
booth #1810
for a demo!



DiaDent Group International
www.diadent.com
1.877.342.3368