

NEW

**DiaDent**<sup>®</sup>

DiaDent Group International  
1.877.342.3368  
www.diadent.com

# DiaRoot<sup>®</sup> BioAggregate

## Root Canal Repair Material



DiaRoot BioAggregate Root Canal Repair Material is a **biocompatible pure white powder** composed of ceramic nano-particles. Upon mixing DiaRoot with BioA Liquid, the hydrophilic BioAggregate Powder promotes cementogenesis and forms a hermetic seal inside the root canal.

### INDICATIONS

- Repair of Root Perforation
- Repair of Root Resorption
- Root End Filling
- Apexification
- Pulp Capping

### FEATURES

- Aluminum-Free Composition
- Biocompatible
- Optimum Working Time & Setting Time
- Pure White Teeth-Colored Powder
- Excellent Workability and Handling Properties

Contact your personal dealer sales representative for  
further information and pricing today!

PRECISION. PURITY. RESULTS