

# DiaDent®

Serving you worldwide



“How many instruments do you want to hold in **YOUR** hand?”



## Miracle Suction

Multifunctional Mirror with Suction Function

A unique and innovative idea!  
These mirrors feature **suction holes** to free up your hands.

- Fewer instruments in the oral cavity
- Less obstruction for the dentists and dental hygienists
- Increased level of comfort for the patients
- Cleaning brush for the hose to remove debris
- 6 holes in the mirror head - rapid suction of liquid while preventing suction on cheek and tongue
- Colour coded handles for easy identification
- Fits into all suction unit on unit chairs
- 100% stainless steel
- Effective especially in pedodontic treatment, scaling and tooth reduction

To place an order, contact your local dealer. For further information, call toll free 1.877.DIA.DENT (1.877.342.3368).