



We'll bend **BACKWARDS**
just for you!

NiTi Endo Hand Files

Increased Flexibility


Five times the flexibility vs. stainless steel files.

Strength

Ten times the stress resistance vs. stainless steel files.

Convenience

Produces smooth and progressive enlargement of root canal. NiTi follows existing canal path and has memory so pre-bending is not necessary.



Efficiency

Color coded handles and pre-inserted silicone stops.

Safety

NiTi resists spontaneous breakage, reduces fracture, and decreases zipping and stripping.

Easy Identification

Symbols on the handle



K-Files



H-Files



Reamers

Available Sizes: .02 taper #8-80 in 21mm and 25mm.

DiaDent[®]

DiaDent Group International
www.diadent.com

Contact your local dealer for further information & pricing.

