

DiaDent[®]

DiaDent Group International Inc.

www.diadent.com

1.877.342.3368

Visit us at the ADA Show Booth #2638

NEW!

Q-Lite[®] LED Curing Light

FEATURES:

- A lightweight durable **aluminum** hand-piece that dissipates heat quickly
- High intensity output of **1100 mW/cm²** that offers fast polymerization
- Wide wavelength **430 - 490 nm** that covers all popular composites
- 4 types of timer setting modes: **5, 10, 15 and 20 sec**
- Convenient preset curing programs that include **two-step** and **soft-start curing**
- Advanced detachable **lithium-ion battery** that does not build memory effect
- **Cordless** operation ensures complete control and freedom of movement



FREE!!

**Dia-400D Light Curing
Protective Glasses**

**Buy 1 Q-Lite,
Get 1 Dia-400D FREE! (Retail Value \$50.00 US)**

Q-Lite Complete Set includes: Hand-piece with battery, light protector, charger base, hand-piece support, power supply (adapter & cord), and operation manual