

DiaDent[®]

www.diadent.com
1.877.342.3368

Presenting the **FUTURE** of teeth whitening technology...

dazzle[™] **NEW** Advanced Teeth Whitening System

BENEFITS:

- Advanced photochemical mediators
- No sensitivity, no pain, and no discomfort
- Superb color lightening
- Repeatable, after appropriate period
- Fast activation: Saves time for dentists and patients
- No need to purchase specialized lamp equipment
- Create and improve patient loyalty

ADVANCED TECHNOLOGY:

Dazzle[™] Powder Mix (Photo Active Booster):

An exclusive formulation that is designed to efficiently absorb light energy coming from various dental curing lights. The Photo Active Booster accelerates the peroxide's activation and shortens the time of bleaching process.

Visual System:

Thanks to Dazzle's visible activation control that is present in the Dazzle Whitening Gel, precise indication lets you know when exactly the whitening gel has achieved its full activation.

INTENSE ORANGE -----> **TONES OF WHITE**

"Cold" Photochemical Activator:

Since an increase in temperature can potentially harm the tooth and lead to tooth sensitivity, Dazzle's COLD type of photochemical activation ensures there is no significant tooth temperature increase, unlike many other whitening products.



Product #: 1010-100 (Single), 1010-000 (Bulk - 3 Singles)

Contents:

- 1 Dazzle[™] Active Shield: Gingival Protection Gel Syringe
Containing 1.2 ml of Metacrylic Polymer Gel
- 1 3ml bottle of 50% Dazzle[®] Hydrogen Peroxide
(Final concentration in active gel: 35%)
- 1 Dazzle[™] Powder Mix: Photo Active Booster
Containing 1.6 g of thickeners, activators, colorant
- 1 Syringe Piston
- 2 Syringe Extrusion Needles
- 1 Spatula
- 1 Instructions for Use
- 1 Treatment Illustration



From dullness to dazzling brilliance...