




DIA-KONDENSOR

Endodontic Condenser

Nickel Titanium



Stainless Steel

FEATURES

- Double ended: Nickel Titanium and Stainless Steel
- ISO sized: compatible with ISO standard files and obturation material
- Durable and lightweight design to reduce hand fatigue
- Ergonomic design with dimpled finger rest
- Color-coded handpiece for fast and easy identification
- Markings on tip to indicate pulp length
- Gradual .02 taper to enhance tactile feedback to prevent root fracture
- Endo stop included on the NiTi end

Size	Item Number
#35/70 (Green)	551-101
#40/80 (Black)	551-102
#50/100 (Yellow)	551-103
#60/120 (Blue)	551-104
Complete Set (4 pcs)	551-191

DOUBLE ENDED
DOUBLE DUTY

