



# DiaDent<sup>®</sup>

DiaDent Group International  
 www.diadent.com | 1.877.342.3368  
 LIKE us: www.facebook.com/diadent




**NEW!**

## Dia-Pro W™ Gutta Percha & Paper Points



- New from DiaDent! Complements **WaveOne<sup>®</sup> Files**
- Patented inspection machine to check diameter accuracy
- One hand slide open package
- Smooth, rounded, and bullet shaped tips
- Readily conforms to shape of the canal
- Optimal apical seal



### GP/PP REFERENCE GUIDE

IF YOU USE:	TRY:
Profile ISO <sup>®</sup> Rotary Files	Dia-ISO GT /Dia-ProISO GT
WaveOne <sup>®</sup> Files	Dia-Pro W (Small, Primary, Large) <b>*NEW*</b>
ProTaper <sup>®</sup> Files	Dia-Pro T (F1, F2, F3, F4, F5)
Twisted Files (TF)	Dia-T
Profile <sup>®</sup> Series 29 <sup>®</sup> Files	Dia-PRO (o4, o6 Taper)
ProSystem GT <sup>®</sup> or ProFile GT <sup>®</sup> Files	Dia-System GTP (#20, #30, #40)
Profile GT <sup>®</sup> Hand Files	Dia-GTP
Shaping Files	Dia-FIT
Oversized for longer/animal canals	Dia-VET (35mm, 45mm, 60mm)
Auxiliary points for master cone technique	Feathered Tip (extra fine taper)
Cordless Obturation and Compaction System	Dia-Gun <sup>™</sup> and Dia-Pen <sup>™</sup>