

DiaDent®



Corporate Office
#704 - 837 West Hastings Street
Vancouver, BC Canada V6C 1B6
[t] 604.629.0231
[tf] 1.877.DIA.DENT (1.877.342.3368)
[f] 604.669.5850
[e] diadent@diadent.com
[w] www.diadent.com

Miracle Suction

Multifunctional Mirrors with Suction Function

A unique and innovative idea!

These mirrors feature suction on the bottom and front sides to free up your hands and save you time.

- Fewer instruments in the oral cavity
- Less obstruction for the dentists and dental hygienists
- Increased level of comfort for the patients
- Cleaning brush for the hose to remove debris
- 6 holes in the mirror head - rapid suction of liquid while preventing suction on cheek and tongue
- Color coded handles for easy identification
- Fits into all suction unit on unit chairs
- 100% stainless steel
- Effective especially in pedodontic treatment, scaling and tooth reduction

Regular Package

Contains:

- 5 autoclavable mirrors with detachable handle
- 2 suction holders
- 2 hoses, 2 adapters, o-rings and cleaning brush

Mini Package

Contains:

- 2 autoclavable mirrors with detachable handle
- 1 suction holder
- 1 hose and 1 adapter



Miracle