



*Flexibility & Strength,
Combined in one...*

NiTi Endodontic Files

Nickel Titanium Hand Files



Increased Flexibility

Ni Ti shaft is made of premium nickel-titanium with five times the flexibility compared to stainless steel files.

Strength

NiTi has ten times the stress resistance compared to stainless steel files.

Convenience

Produces smooth and progressive enlargement of root canal.
NiTi follows existing canal path and has memory so there is no need to pre-bend the files.

Efficiency

Color coded handles and pre-inserted silicone stops.

Safety

NiTi resists spontaneous breakage.
Increased angular rotation reduces fracture.
NiTi decreases zipping and stripping.

Easy Identification

Symbol on the handle:



K-Files



H-Files



Reamers

DiaDent Group International

#11 - 3871 North Fraser Way
Burnaby, BC Canada V5J 5G6
t: 604-451-8851
f: 604-451-8865
tf: 1-877-342-3368
www.facebook.com/diadent

Available Sizes:

- K-Files (6/box): #8 - #80 in 21, 25 and 31mm
- H-Files (6/box): #8 - #80 in 21, 25 and 31mm
- Reamers (6/box): #8 - #80 in 21, 25 and 31mm



Serving you worldwide.