

NEW

e-Temp™

Hydraulic Temporary Restorative Material



Order through your local dental dealer | Item #2004-1001 (30g white)


Indications

- ★ Temporary restorations - especially inlay restorations
- ★ Temporary filling and sealing
- ★ Temporary filling of implant sites


Features

- ★ Adheres to dentin, expands slightly during setting and hardens rapidly to ensure hermetical seal
- ★ Improved mechanical properties for quick temporary filling
- ★ Moisture-activated: sets within 25 minutes or less
- ★ Contains zinc oxide and calcium sulfate
- ★ Excellent adaptability
- ★ Easy to remove

Setting Time



Compressive Strength



DiaDent®

www.diadent.com | 1.877.342.3368

Scan QR code, LIKE us on FACEBOOK and receive an e-Temp sample!*



*Canadian customers only. Certain conditions apply. See our FACEBOOK page for details.